



“Indego is a therapy tool that allows me to gait train my patients like never before. It’s powerful and easy to use.”

—Casey Kandilakis, PT, DPT, NCS, Shepherd Center, Atlanta, Georgia



Assist-as-Needed Support

Individualized therapy thanks to variable assist functionality. Clinicians can challenge their patients, make adjustments to suit their patient's individual needs, and maximize patient's recovery potential. Various gait parameters can be adjusted in real-time to enable a physiological gait pattern and goal-based motor skill training. Indego offers clinicians a gait therapy tool which allows for unprecedented training variability on various surfaces.

Clinical Applications Expertise. Cutting edge technology in combination with outstanding clinical expertise is a recipe for success. That's why we focus not only on the development of superior technology, but also on clinical excellence. A comprehensive Indego training concept with various competency levels will allow therapists to become Indego experts and to leverage the enormous therapy potential offered by Indego.



“ High patient satisfaction, a modern therapy offering and enthusiastic and efficient therapists are of highest importance to me – Indego helps to accomplish this. ”

—Clair Jones, MS, Barrow Neurological Institute, Rehabilitation Services, Phoenix, Arizona

Fast Setup to Maximize Efficiency

Focus on walking and not on device setup. Thanks to its modular design, a fast patient setup is guaranteed. Therapists and patients appreciate a quick learning curve thanks to Indego's intuitive approach. Therefore, valuable training time is spent on gait therapy and not on device learning and setup.

Three devices - treat more patients. The Indego Therapy Kit comes in three devices, which increases productivity in your clinic. More steps per session can be achieved with less physical exertion on the part of the therapist compared to traditional gait training. In addition to improving staff efficiency, this cutting edge technology can attract patients and referrals. Healthcare administrators are also attracted by a fast return on investment.


“Indego is an extremely modular exoskeleton which is very easy to use, thus giving patients the ability to walk with it during the first session.”

—Arun Jayaraman PT, PhD,
Rehabilitation Institute of Chicago,
Chicago, Illinois

Based on Motor Learning Principles

From clinicians for clinicians. Functional, task-specific, over-ground training in combination with high patient engagement and the provision of performance feedback enables therapists to offer gait therapy which is based on scientifically proven motor learning principles.

Parker is committed to developing future Indego features derived from evidence based findings. Working closely together with experts in the field, Parker strives for excellence and superiority.



Data from the Indego multi-center study. Figure shows average results of 45 subjects from the 10 Meter Walk Test - Outside. Subjects were tested midway through and at the end of the study. Significant improvements in walking speed have been shown. Highest walking speed measured was 0.59 m/s among all subjects. Subjects continued to improve gait speed through the end of the study even after becoming proficient with the Indego midway through the study.







indegocom



Shepherd Center

Shepherd Center is Parker's lead clinical partner for Indego.

FDA Clearance
for Clinical and Personal Use



Parker Hannifin Corporation, Human Motion & Control,
1390 E. Highland Rd. Macedonia, OH U.S.A.
1.844.846.3446 support.indego@parker.com

The Indego orthotically fits to the lower limbs and the trunk; the device is intended to enable individuals with spinal cord injury at levels T7 to L5 to perform ambulatory functions with supervision of a specifically trained companion in accordance with the user assessment and training certification program. The device is also intended to enable individuals with spinal cord injury at levels T4 to T6 to perform ambulatory functions in rehabilitation institutions in accordance with the user assessment and training certification program. The Indego is not intended for sports or stair climbing.



Indego® is a registered trademark of Parker Hannifin Corporation. © 2017 Parker Hannifin Corporation