

evolvrehab[®] CE

Play. Enjoy. Improve!



evolv



EvolvRehab® is a CE certified **virtual rehabilitation platform** which has been used around the world to provide **clinical** and **telerehabilitation services**.

Patients can perform **therapeutic activities** in clinics and the comfort of their own homes through **personalized exergames**, ROM assessments and a variety of exercises routines.



Therapists use the EvolvRehab® Manager to **remotely set up** patient profiles, **prescribe treatment plans** consisting of personalized virtual therapy activities, **monitor patient adherence and results** and communicate with patient via in-app messages or through video conferencing.



EvolvRehab® can be used to help treat a wide **variety of neurological and age-related conditions and diseases** including stroke, MS, Parkinson's, Cerebral Palsy, TBIs and frailty in older adults.






Photo courtesy of Microsoft

The RehabKit only needs to be connected to a power outlet and a TV to begin providing EvolvRehab[®] therapy in clinics and at home. The **simple one-button start system** means **no technology experience is required**.


The Evolv RehabKit **includes all the hardware necessary** to provide virtual rehabilitation services in clinics and in patients' homes. It incorporates the **Azure Kinect 3D camera** for body tracking and a powerful **gaming mini pc**.



The RehabKit allows patients to have **greater autonomy** when choosing when to do their rehabilitation, while allowing them to maintain close contact with their therapists to **ensure adherence**.

Advantages of EvolvRehab® For Clinics & Telerehabilitation

For Therapists



- ✓ Setting up of patient profiles
- ✓ Personalization of virtual therapy sessions
- ✓ Tracking of patient adherence + performance
- ✓ Video consultation check-ups via Teams



The Evolv virtual rehabilitation package includes the Evolv RehabKit + the EvolvRehab® On Premises or Cloud Station Subscription License

For Patients



- ✓ Access to home-based therapy
- ✓ Goal setting
- ✓ Simple to use hardware
- ✓ Available 24-7 365 days/year
- ✓ Video consultations via Teams

How EvolvRehab Works in Hospitals, Clinics, Residential Care and Home Settings

