

**LANDICE**



# SPORTS/PHYSICAL THERAPY & REHABILITATION

**REHABILITATION TREADMILLS - ELLIPTICALS - BIKES**

The Landice line of rehabilitation products provides the very best in functional usage and reliability. From reversing treadmills to low-impact ellipticals and easy-to-start self-powered bikes, Landice fitness equipment provides the perfect tools to assist a wide variety of patients in their recovery and rehabilitation.

## **REHABILITATION EQUIPMENT**

- L7 / L8 Rehab Treadmills
- E9 Elliptical
- R9 / U9 Recumbent & Upright Bikes



# THE ROAD TO RECOVERY STARTS HERE

## L7 & L8 REHAB TREADMILLS

Made with the finest materials, the L7 and L8 Rehabilitation Treadmills provide the best and safest support for all your rehabilitative programs. With a low step-up height, zero starting speed, 300 micro-amp isolation leakage kit and extended parallel medical handrails, these treadmills are ideal for any physical therapy and sports conditioning center. A number of unique options are also available, such as a reversing treadbelt for downhill walking, 100 micro-amp hospital isolation leakage kit, and an orthopedic suspension system. For a complete list of options please visit us at [Landice.com](http://Landice.com).

## TREADMILL FRAME COMPARISON

Solid aluminum construction - rust-free and durable

	L8	L7
<b>User Capacity</b>	500 lbs.	400 lbs.
<b>Equipment Weight</b>	355 lbs.	315 lbs.
<b>Dimensions</b>	36" W x 83" L x 61.5" H	36" W x 76.5" L x 60" H
<b>Running Surface</b>	22" W x 63" L	20" W x 58" L
<b>Power</b>	110 VAC, 60Hz, 15A 220 VAC, 60Hz, 15A (optional)	110 VAC, 60Hz, 15A 220 VAC, 60Hz, 15A (optional)
<b>Rollers</b>	3 1/2" 22 lbs.	2 1/2" 14 lbs.
<b>Step-up Height</b>	7 1/2"	5 1/2"
<b>Speed</b>	0.1 to 12 mph	0.1 to 12 mph
<b>Incline</b>	0% to 15%	0% to 15%
<b>Drive System</b>	4 HP Continuous Duty Motor	4 HP Continuous Duty Motor
<b>Deck</b>	1" - Reversible	1" - Reversible
<b>Control Panel</b>	Rehabilitation	Rehabilitation
<b>Heart Rate</b>	Contact Heart Rate	Contact Heart Rate
<b>Warranty*</b>	5 yrs. parts/2 yrs. wear items/ 1 yr. labor	5 yrs. parts/2 yrs. wear items/ 1 yr. labor

\*Treadmills are covered by the LTD warranty which offers coverage for usage of less than 5 hours per day.



## COMFORTABLE AND SAFE REHAB

### BIKE SPECIFICATIONS

	R9	U9
Frame Type	Step-through recumbent with low step-over height	Upright
Pedal Type	Self-balanced pedals with adjustable straps	Self-balanced pedals with adjustable straps
Seat Adjustment	Forward/backwards with 15 settings	Up/down with 12 settings
Flywheel Weight	15 lbs.	15 lbs.
Resistance System	Eddy current magnetic braking system	Eddy current magnetic braking system
Resistance Range	1-20	1-20
Handlebar	Console grip bars, seat grips with contact heart-rate monitoring and resistance controls	Racing style handlebars with contact heart-rate and resistance controls
Display	LED with multiple feedback windows	LED with multiple feedback windows
Heart Rate	Contact heart rate for heart rate control/ wireless chest-strap (optional)	Contact heart rate for heart rate control/ wireless chest-strap (optional)
Accessories	Water bottle holder, reading rack, accessory tray	Water bottle holder, reading rack, accessory tray
Cooling Fan	3 speeds	3 speeds
Power Requirement	Self-powered	Self-powered
Transport Wheels	Yes	Yes
Levelers	Yes	Yes
Assembled Weight	205 lbs.	159 lbs.
Dimensions	65.74" L x 25.31" W x 52.20" H	44.56" L x 24.72" W x 58.11" H
Max User Weight	350 lbs.	350 lbs.
Warranty*	5 yrs. parts/1 yr. labor/90 days for high-wear items	5 yrs. parts/1 yr. labor/90 days for high-wear items

\*Bikes are covered by the Club warranty which offers coverage for unlimited hours of usage.

### R9 / U9 BIKES

The Landice R9 & U9 bikes offer a comfortable and safe rehabilitation option. The wide step-through design on the R9 allows for easy access along with superior back support. The U9 offers well-cushioned seating with adjustable height levels to support a personalized lower-body rehabilitation regimen.



## NATURAL BIOMECHANICS - LOW IMPACT

Patients at any athletic level from beginner to professional will find physical benefits and rehabilitative functionality on a Landice.

The E9 offers one of the only fitness units that enables both upper and lower body exercise capabilities. Patients can stand upright on the frame's stationary footpads and perform an isolated workout to strengthen the upper body. This back-and-forth motion is considered a preferable position as opposed to a seated rotational movement.

For the lower body, the center-drive technology replicates the biomechanics of walking or running. The E9 provides a low impact alternative to treadmills.

The elliptical also features an accessible, safe, stable and patient-friendly exercise platform. Users can easily grab the side rails, step-up and safely start a workout. The pedals will move only when the patient engages the upper body handrails and/or foot pedals.



### E9 ELLIPTICAL FEATURES

<b>User Capacity</b>	500 lbs.	<b>Stride Length</b>	Fixed 21"
<b>Active Footprint</b>	35" W x 76" L x 67" H	<b>Drive System</b>	Center Drive
<b>Physical Footprint</b>	35" W x 68" L x 67" H	<b>Flywheel</b>	11" diameter, 20 lbs.
<b>Power</b>	110 VAC, 60Hz, 15A 220 VAC, 60Hz, 15A (optional)	<b>Panel</b>	Cardio
<b>Handrails</b>	2" Medical-Grade Side Rails	<b>Pedals</b>	2" Spacing with Orthopedic Gel Inserts
<b>Warranty*</b>	5 years parts/2 years wear items/1 year labor		

**Certifications: UL 1647, CAN/CSA - C22.2, CE**

\*Elliptical is covered by the Club warranty which offers coverage for unlimited hours of usage.

For more information on E9 Cardio, please visit us at [Landice.com](http://Landice.com).

Landice engineers, assembles and tests the finest specialty fitness equipment on the market and is one of the only manufacturers to do so in the USA. Landice has been recognized repeatedly by leading consumer rating publications as having the world's finest treadmills and ellipticals.

Visit us at [www.landice.com](http://www.landice.com) for additional options and complete warranty information.

