

**TOPRO
TAURUS**

TOPRO



STEADY!

Safe walkers for institutions and at home

www.topro.co.uk

TOPRO TAURUS

Steady and safe walkers for institutions and at home

- Easy to read heightscale and easy to adjust the height.
- High enough for tall users (up to 2 m).
- Tested and approved according to ISO 11199-3 for use indoors on plain foundation.
- The TOPRO TAURUS E-models are tested and approved according to IEC 60601-1:2005 + corr 1:2006 + corr 2:2007 and IEC 60601-1-2:2007.

Maintenance free material

- Aluminium with polyester surface and polyurethane.
- Easy to clean.
- Non-corrosive.
- All wear and tear parts can be exchanged.

Two height adjustment methods:


TOPRO TAURUS H Hydraulic height adjustment.

TOPRO TAURUS E Electrical height adjustment.

Two levels of equipment:

BASIC Standard equipment with adjustable arm pads and foot brakes on all wheels. Delivered partly assembled.

PREMIUM Premium has in addition:



	Art. No.
TOPRO TAURUS H BASIC	814785
TOPRO TAURUS H PREMIUM	814787
TOPRO TAURUS E BASIC	814790
TOPRO TAURUS E PREMIUM	814789

TOPRO TAURUS H

- Gas springs simplify the height adjustment.
- Can be used in the shower.

Product weight 18 kg.

Maximum weight of user 150 kg.





TOPRO TAURUS E

- Electric motor makes the height adjustment very simple.
- Can also be used as raising support.

Product weight 22 kg.

Maximum weight of user 150 kg.



Foot plate

814645



Forearm side support (pair)
(Launch winter 2014/15)

814647



Brake handle
left side

814846



Adjustable handles

814001



Tray

814810



Height adjustment panel
(Concens)

814646



Wall mount for the
Taurus battery charger (Concens)

814085



Small basket

814821



Infusion stand

814793



Stand for
body fluid-boxes

814034



Oxygen bottle
holder

814820