

# Levitop Comfort eFlow & Standard eFlow

**NEW**

Mobile Patient Lift

Origin: Poland  
Certificate: CE



**Introducing Levitop Comfort eFlow and Standard eFlow, a patient lift that improves daily patient care. This solution increases comfort and safety for both patient and caregiver.**

**The Levitop eFlow is equipped with a system that assists in the movement of the lift. This allows the operator to easily guide, rotate and turn the lift with the patient using minimal effort.**

## Features

- The lift is equipped with a four-point sling as standard
- The sling rotates 360° with the possibility of locking in four positions (90°)
- Wheels allowing maneuvering even in small and narrow spaces
- Rear wheels equipped with a brake
- Remote control included
- Adjustment of the width of the leg span with the remote control
- Emergency lowering mechanism to increase patient safety

## Levitop Comfort



## Levitop Standard



## eFlow Drive System



### Drive System

The drive system, powered by LINAK's lithium-ion battery, allows the lift to be moved in any direction with minimal effort



### Steering Handle

The steering handle mimics the operator's movements thanks to built-in sensors

## Specifications

	Standard eFlow	Comfort eFlow
Frame length (cm)	150	150
Frame width (cm)	67	67
Device height (cm)	135	135
Push bar height (cm)	92.8	92.8
Raising height (cm) (min./max.)	68/196	68/196
Minimum legs width (cm) (int./ext.)	52/67	52/67
Maximum legs width (cm) (int./ext.)	103/118	103/118
Turning diameter of the product (cm) (min)	144	144
Turning diameter of the product (cm) (max)	158	158
Lifting speed (s) (raising)	36	36
Lifting speed (s) (lowering)	34	34
Diameter of castors with brake (cm)	12.5	12.5
Diameter of castors without brake (cm)	10	10
Safe working load (kg) (SWL)	≤230	≤230
Device weight (kg)	61	66

# VertiUp eFlow & VertiUp Plus eFlow

Interactive Wall Projection

**NEW**

**Origin:** Poland  
**Certificate:** CE Class I



**VertiUp / VertiUp Plus eFlow patient lift is designed to assist in transferring and repositioning individuals with reduced mobility due to illness or disability. It has a wide range of applications in assisting personnel in caring for patients affected by conditions such as craniocerebral trauma, spinal cord injury, disease, and Parkinson's disease, as well as multiple sclerosis.**

### Features

- Can be used in small spaces
- Emergency lowering function
- Maximum safety and comfort
- Overload protection and very easy to maneuver
- Relieves the strain on the operator's spine during handling
- Ability to drive under a chair, wheelchair or toilet, thanks to the device's swinging legs
- Electric control of the remote control for raising, lowering and spreading the legs

### VertiUp



### VertiUp Plus



## eFlow Drive System



### Drive System

The drive system, powered by LINAK's lithium-ion battery, allows the lift to be moved in any direction with minimal effort



### Steering Handle

The steering handle mimics the operator's movements thanks to built-in sensors


## Specifications

	VertiUp	VertiUp Plus
<b>Base legs spreading:</b>	electrical	electrical
<b>Max. safe working load (kg) (SWL):</b>	250	310
<b>L x W x H (cm):</b>	135 x 68.2 x 188 - 225	143.3 x 77.9 x 188 - 225
<b>Lifting range (cm):</b>	34.1 - 172.6	32.2 - 170.6
<b>Min. turning diameter (cm):</b>	170	184
<b>Min. base legs width (cm) (int./ext.):</b>	53.5 / 68.2	63.7 / 77.9
<b>Max. base legs width (cm):</b>	125.9 / 140.6	143.3 / 152.8
<b>Base legs height (cm):</b>	11.1	7.3
<b>Wheels diameter (cm):</b>	10 / 12.5	5 / 12.5
<b>Weight (kg):</b>	87.7	97.6
<b>Battery capacity</b>	2.9 Ah	2.9 Ah
<b>Remote control function:</b>	<ul style="list-style-type: none"> <li>• up / down</li> <li>• hanger bar forward / backward</li> <li>• el. leg spread: outside / inside</li> </ul>	<ul style="list-style-type: none"> <li>• up / down</li> <li>• hanger bar forward / backward</li> <li>• el. leg spread: outside / inside</li> </ul>

**Origin:** China

A dark grey leather recliner chair with a footrest extended, set against a blue background. The chair features a high back, padded armrests, and a large storage pocket on the side. The footrest is extended forward, and the chair is mounted on a black metal base with casters.

Height (cm)	117
Depth (cm)	90
Width (cm)	75
Seat Height (cm)	48
Seat Depth (cm)	52
Seat Width (cm)	50
Back Cushion Height (cm)	76
Max User Weight (kg)	130
Li-on Battery (mAh)	10000



# Ravenna

**NEW**
**Origin:** China

- The seat can tilt forward to help elderly users safely get off the chair.
- The headrest and lumbar cushion feature electric air inflation for optimal support.
- Detachable armrests; optional detachable dining table available to increase flexibility and caregiving convenience.
- Extended footrest and reclining backrest design enhance overall comfort.
- Equipped with manual wheels and a backrest push handle for easy maneuverability.
- IV pole can be added to meet medical care needs.



## Specifications

Unit Height (cm)	104
Unit Width (cm)	76
Unit Depth (cm)	89
Unit Weight (kg)	55.75
Seat Height (cm)	47
Seat Depth (cm)	51
Seat Width (cm)	40/51
Back Cushion Height	65
Max User Weight (kg)	130


**Full Relax**

**Lift**

## Vision Electric Standing Frame

Standing Frame With Power Lift

Origin: Germany  
Certificate: CE, ISO 60601-1

### Features

- Wireless remote control
- Easy to use height adjustment of the therapy table by means of gas springs
- Depth and tilt adjustable therapy table
- Adjustable depth, height, breadth and angular position of knee supports
- Width and depth adjustable heel supports
- 4 safety wheels



### Optional Accessories


- Various belt systems
- Posterior support
- Back support with headrest
- Foot fixation, Velcro
- Forefoot fixation

### Specifications

Standing Frame		Vision® Junior	Vision®
Article Number		VSJ-1000	VS-1100
Footplate Table	cm	87 to 100	93 to 120
Footplate Knee Support	cm	30 to 50	38 to 56
Height Of Handrails (From Footplate)	cm	69 to 90	80 to 104
Floor Upper Edge Footplate	cm	8	8
Chassis Length	cm	93	103
Overall Length	cm	100	110
Outer Frame Width	cm	68	78
Inner Frame Width	cm	51.5	61.5
Therapy Table (W x D)	cm	65 x 60	75 x 60
Depth Adjustment Table	cm	16	16
Tilt Adjustment Table	Degree	+/- 15°	+/- 15°
Unit Weight	kg	58	62
Maximum Load	kg	150	150



Raise up using a pressure relieving stand up belt




Adjustments at the standing frame with power lift Vision®



3-D adjustable pelvic support (accessory)



3-D adjustable knee supports



Standing frame with power lift Vision® and wireless remote control



Quality with safety



## Prelivia **NEW**

Neurostimulation Device for Preventing Pressure Injuries


**Origin:** Canada  
**Certificate:** FDA

**Prelivia is aimed at protecting patients from pressure injuries, also known as pressure ulcers and bedsores. It uses proprietary neurostimulation that is scientifically proven to continuously stimulate blood circulation and minimize tissue damage.**



	Prelivia
Temporary pressure relief	✓
Pressure redistribution	✓
Restore blood circulation	✓
Increase tissue oxygenation	✓

### Prelivia is Scientifically Proven



1 – Gyawali, et al., J Apple Physiol. 110:246-255, 2011;

2 – Solis et al., J Appl Physiol, 114: 286-296, 2013;

3 – Ahmetovic, et al., Adv Wound Caref, 4(3): 192-201, 2015. Kane, et al, Adv Wound Caref, 6(4):115-124,2017;

4 – Safety and Acceptability Study, Chan, et al., 2017



## Supro Tilt Table

A Unique Supine Stander  
That Has Been Designed For Children From 2 Years To Young Adult

Origin: UK

The Supro provides the option of supine, prone and vertical standing in one frame. Its design gives an angle range of 90° allowing it to achieve a horizontal and vertical position.



### Features

- Supine, prone and vertical standing in one
- Highly adjustable
- Supportive and comfortable positioning
- Support to the knee, hip and chest
- Politeal pads provide the necessary comfort
- Range of accessories available
- Available in choice of 2 widths and 2 sizes

### Supro Accessories


The Supro comes complete with upholstered board, lateral and pelvic supports, knee cups, cocoon sandals and plain tray.



Pelvic/Trunk Lateral Supports



Soft Head Support



Padded Tray Insert



**Kneeling Box**  
Enables use in an upright kneeling position.



**Knee Cups**  
Maintains knee in position  
Available in sizes 1 - 6



**Sheepskin Liners For Kneecaps**  
Maintains position of knee with added protection  
Available in sizes 1-3



**Popliteal Foams**  
Provides support to back of knee



**Knee Cap Spacers**  
Provides additional height to the popliteal pad to accommodate differing degrees of knee contraction

## Heathfield Chair

Fits Well Into Nursery, Clinic, School & Home Environments

Origin: UK

### Features

- Perfect starter chair
- Adjustable seat height and depth for growth
- Ideal for nursery, school and home
- Extremely durable
- Wide range of accessories
- Available in 6 sizes



### Heathfield technical data all dimensions in mm unless otherwise stated

Adjustable Arm Code	Fixed Arm Code	Size	Seat Height*	Width	Seat Depth	Weight Limit (kg)	Back Height
6901	6911	1	180 - 230	240	180 - 240	25	250
6902	6912	2	200 - 300	280	220 - 310	45	280
6903	6913	3	240 - 370	325	260 - 340	50	320
6904	6914	4	290 - 390	365	320 - 400	60	365
6905	6915	5	325 - 425	435	370 - 460	70	420
6906	6916	6	395 - 545	435	450 - 550	85	420

### Accessories



#### Skis

Provides greater stability



#### 4 Glisdomes

Comes complete with screws to fix under chair leg



#### Reversible Footboard

Foot support that prevents unwanted chair movement



#### Footplate & Straps

Ensures feet are firmly in position



#### Kicking Board

Maintains foot position  
Strongly recommended for a child who straps leg under the chair



#### Mobile Dolly Base

Facilitates easy manoeuvrability



#### Adjustable Height Push Handle

Adjusts to achieve optimum push position. To be used in conjunction with a mobile dolly base or minimum raise mobile base



#### Pelvic Cushions - pair

Provides hip and thigh support  
Available in various sizes



#### Pommel

Encourages hip abduction



#### Adjustable Footrest

The sides help to keep feet contained and the footrest is easy to adjust



#### Tray Handrail Full Width

Helps to stabilise upper limbs  
Additionally can be used to attach toys



#### Wooden Tray

Promotes trunk support  
Additionally can be used for activities, playing and feeding

## POPULAS Therapy Table

Origin: USA

### 1 Adjustable Worktable With Tilt

- Adjustment range 23" -33"
- Tilt adjust from 0° to 50°
- Oversized knob adjusts in 1-inch increment

Product Code	Dimension (L X W)
KA-3630-T	36" X 30"
KA-3630	36" X 30"



### 2 Accella Adjustable Workstation/Table

- Adjustment range is 27" -39"
- Hand-crank adjusts at 5 turns per inch
- Laminate top 1-1/8" thick
- Heavy 3mm PVC edge banding

Product Code	Dimension (L X W)
AC 7236	72" X 36"
AC 4836GT	48" X 36"
AC 6648	66" X 48"



### Accessories



ADA Signage



Casters



Locking Casters



Comfort Curve

Cable Management

### 3 Equity Tilting Workstation

- Adjustment range 24" -38"
- Hand-crank adjusts at 6 turns per inch
- Durable laminate top or real wood veneer
- Handle retracts under table out of user's way
- Tilts from 0-60 degree

Product Code	Dimension (L X W)
EQ 3630CR-T	36" X 30"
EQ 4830CR-T	48" X 30"



# Changing & Showering Benches

Wall Mounted & Space Saving For Dry Or Wet-Room Environments

Origin: UK

## 1 Hi-Riser Shower Stretcher

- Both changing and showering
- Wall mount and floor mount
- Loading up to 150kg
- Electric height adjustable
- Stainless steel construction
- Water proof remote control
- Excellent height range for easier handling
- Quick drainage
- Available in 2 sizes



Stainless Steel	Size	Table Dimensions			Height Range	Wall Frame Width	Wall Frame Depth	Wall Frame Height	Safety Side Height	Stored Width	User Weight	Shipping Weight
		L mm	W mm	H mm								
3063	1	1400	650	1110	390 - 960	700	110	1550	165	300	150	80
3067	2	1800	650	1110	390 - 960	700	110	1550	165	300	150	85

## 2 Hi-Riser Changing Bench

- Wall and floor mounted
- Excellent height range
- Stowaway side rail
- Easy to clean padded surface
- Choice of colours and sizes
- Frames available for both dry and wet-room environments
- Wall mounting kit included



Mild Steel Powder Coated	Stainless Steel	Size	User Weight  kg	Table Dimensions			Travel Range  H mm	Wall Frame Width  A mm	Wall Frame Depth  B mm	Wall Frame Height  C mm	Safety Side Height  Y mm	Stored Width  S mm	Max User Weight  kg	Product Weight  kg
				L mm	W mm	M mm								
3041	3051	1	200	1200	550	875	400 - 1000	700	110	1550	165	305	200	80
3042	3052	2	200	1300	550	875	400 - 1000	700	110	1550	165	305	200	82
3043	3053	3	200	1400	600	925	400 - 1000	700	110	1550	165	305	200	84
3044	3054	4	200	1500	600	925	400 - 1000	700	110	1550	165	305	200	86
3045	3055	5	200	1600	650	975	400 - 1000	700	110	1550	165	305	200	88
3046	3056	6	200	1700	650	975	400 - 1000	700	110	1550	165	305	200	90
3047	3057	7	200	1800	650	975	400 - 1000	700	110	1550	165	305	200	92



# TOPRO Taurus E Premium Walker

Stable Walker For Home Or Care Facility

Origin: Denmark  
Certificate: CE, IEC 60601-1

## Features

- Fully assembled with adjustable handles, driving and parking brake on the righthand handle and the opportunity to lock the swing (directional block)
- Easily visible red button for height adjustment which can be connected to external height adjustment panel (supplementary equipment)
- Adjustable distance between forearm supports
- Easy to read height scale and an electric motor to make the height adjustment very user friendly
- Central column provides easy access for nursing staff
- Foot brakes on all four wheels
- Driving and parking brakes can easily be used by the righthand handle, and work on both rear wheels simultaneously
- Handles can be adjusted forward and sideways for best fit for the user
- Locking the swing (directional block) is possible for the rear wheels



## Specifications

Total Length	78cm
Total Width	66cm
Handles	Adjustable
Forearm Supports	
Length	34cm
Width	34cm to 47cm
Height	86.5cm to 136cm
Diameter	
Handles	3.5cm
Turning Circle	91cm
Wheels	
Diameter	12.5cm
Width	3.1cm
Brake	
Driving and parking brake on the handles	Yes
Foot Brakes	Yes (All wheels)
Directional Block	Yes
Unit Weight	22kg
Maximum Load	150kg
Recommended User Height	140 - 210cm
Product Code	814789

## Accessories



Foot Plate



Stand For Body Fluid-Boxes



Tray



Height Adjustment Panel

## Eleveo

System For Dynamic Unloading

Origin: Poland  
Certificate: CE



**Eleveo is a training device for simultaneous control of unloading, posture and balance on a treadmill or firm surfaces. Systems providing dynamic unloading are a perfect solution for training patients in a wide range of gait disorders.**

**Eleveo helps to maintain proper posture, reduces load, eliminates balance problems and improves motor coordination training. Unique harnesses are designed for uni- or bilateral support of the body as well as for a gradual level change from fully loaded to fully unloaded. The device allows for a manual controlled treatment of lower limbs and pelvis to ensure proper movement. Configurability of the device allows the clinics to treat various patients with a single Eleveo system. Extra options help to adjust the system to specific needs of each patient.**



Swivel sling  
(for sideways gait)



Electronic  
measurement unit



## Assistive Aids

### Features

- Battery-powered electric height adjustment (charged from a 230v wall outlet)
- Two-point suspension with pelvis positioning and front-back inclination using 4 belts
- Electronic measurement unit for monitoring lateral and overall unloading level and visual feedback for the patient
- Four-wheel system with two wheels with full and two with directional brake (the wheels with directional brake allow for pre-setting the direction of movement prior to the treatment)
- The system can lower to 164cm for an easy access to rooms with low doors (180cm) or treatment of shorter patients
- 89.4 cm wide frame allows easy passage through most doors
- Maximum height of 234cm
- Patient height of up to 210cm
- Patient weight of up to 160kg
- Front, back and sideways gait training
- Gait direction change possible without disconnecting the harness
- Variable angle adjustable handlebars
- Dynamic unloading system which allows for shifting the centre of gravity by at least
- 5cm for a more physiological movement of the body during the training
- Including 2 universal suspension harnesses M (size: 92-100 cm), XL (size:132-142 cm)  
(made from washable fabric for easy cleaning)



### Specifications

<b>Dimensions (L x W x H)</b>	135cm x 89.4-109.4cm x 164-234cm
<b>Unit Weight</b>	120kg
<b>Maximum Load</b>	160kg

\*The gate is factory-set to dimension of 730mm.

\*The gate construction allows it to be set up from 730 to 930mm with 25mm movement- settings in total.



### Related Items

- Used for training with Eleveo or stationary hoists

## Neater Eater Robotic

Enables People To Feed Themselves

Origin: UK  
Certificate: CE

### Features

- Controlled by a touchscreen or plug-in switches
- Touchscreen menus and help guide you through many options such as setting mouth position, automatic spoon-wiping and plate turning
- Removeable Touchscreen controller with on-screen help and video guides
- Robot Arm has wide range of movement to reach user's mouth



## Neater Eater Manual

Enables People To Feed Themselves

Origin: UK  
Certificate: CE

### Features

- Different attachments and settings are suitable for people with different abilities
- Available both right and left handed
- Plate flick guard option stops spoon from spilling food at edge of plate
- Plate holder fitted in different positions for accessibility
- Clamps standard, compact



### Benefits

<b>Functional</b>	<ul style="list-style-type: none"> <li>Enable users to eat, drink, call, scratch their head, wipe their nose on their own</li> <li>Supports the user's arms when using keyboard</li> </ul>
<b>Physical</b>	<ul style="list-style-type: none"> <li>Helps maintain strength</li> <li>Reduces development of postural deformities</li> <li>Increases movement range, reduces the chance of developing joint contractures</li> </ul>
<b>Psychological</b>	<ul style="list-style-type: none"> <li>Depends less on carers</li> <li>Encourages social interaction</li> <li>Helps maintain dignity</li> <li>Promotes engagement in activities and mealtimes</li> </ul>

### Saebo Glove

Independence At Your Fingertips

Origin: USA  
Certificate: CE

The proprietary tension system extends the fingers and thumb following grasping so clients can incorporate their affected limb functionally. This new found freedom leads to improved motor recovery and functional independence.



#### Features

- Spiral forearm design that secures the wrist in a functional position
- Tensioners are located at the IP joint (interphalangeal joint) of the fingers and thumb to assist with extending the digits following grasping
- Individual tensioners can be removed to customize assistance based on the client's needs
- Numerous sized tensioners included to accommodate various finger lengths
- Full joint finger motion possible to maximize functional performance
- Non-slip liner to minimize migration



### SaeboStep

Walk Smarter. Confidence And Comfort Is One Step Away

Origin: USA  
Certificate: CE

SaeboStep consists of lightweight, uniquely designed foot drop brace that provides convenience and comfort while offering optimum foot clearance and support during walking.

#### Features

- **Lift:** Strong durable spectra cord easily slips onto included eyelet attachments and lifts the foot quickly and easily
- **Adjust:** The revolutionary BOA dial technology allows individuals to quickly customize the lift angle required for safe foot clearance during walking
- **Secure:** Hook and loop Velcro strapping system secures the device to the ankle
- **Release:** Conveniently release tension at any given time during the day with the simple-to-use dial technology
- **No Laces. No problem:** With the accessory kit to enable shoes without eyelets to be modified, the SaeboStep can even be worn comfortably with the majority of male or female shoe styles



#### Updates

- A new ankle cuff dial and magnetic fastener for quick and easy one-handed donning
- Sleek and sturdy design for added stability and long-lasting comfort
- Redesigned shoe eyelet hooks for less wear and tear and improved performance

## SaeboMAS

Maximize Function While Improving Strength & Control

Origin: USA  
Certificate: CE

**The SaeboMAS dynamic mobile arm support system is a zero gravity upper extremity device specifically designed to facilitate and challenge the weakened shoulder and elbow during functional tasks and exercise drills.**

### Features

- Adjustable spring based parallelogram offering various levels of assistance
- Incorporates all planes of movement so multi-directional activities can be achieved
- Measurable graded tension scale for tracking and documenting progress
- Height adjustable
- Table mount
- Comfortable malleable forearm support with removable liners for infection control
- Elbow support attachment
- Lightweight and portable



## SaeboMAS mini

Unweight The Arm & Explore Life

Origin: USA  
Certificate: CE

**The SaeboMAS mini is the ideal home program solution to assist with improving strength, motor recovery, and independence. Clients will be able to perform functional exercises with greater ease and minimal compensation. Additionally, this personal device allows for enhanced independence for self-care, leisure or occupational tasks like using a computer, eating, drinking, or grooming.**

### Features

- Adjustable spring based parallelogram offering various levels of assistance
- Incorporates all planes of movement so multidirectional activities can be achieved
- Height adjustable
- Includes a portable table mount
- Comfortable malleable forearm support with removable liners for cleaning
- Lightweight and portable





### Carbonhand 2.0 **NEW**

The Smart and Intuitive Grip-strengthening Glove

Origin: Sweden  
Certificate: CE

**Carbonhand is an assistive device with pressure sensors that detect and amplify the user's grip. The more force the user applies, the more Carbonhand enhances it. Designed for daily activities, it can be used indoors, outdoors, and in public places. Worn like a regular glove, it connects to a compact power unit worn close to the body.**



#### Features

- A three-finger glove that enhances the grip strength during usage
- Multiple sensors on each finger and one on the first MCP joint
- Control Pad with three program- mable profiles and an activation button, configured via smart app
- Various glove sizes and carrying options to fit individual needs
- Proven to enhance the user's ability to perform activities in daily living.



#### Enhances Grip and Endurance

Artificial tendons and electric motors dynamically add power to the glove ensuring a firm grip and endurance. It supports daily activities, work, or rehabilitation for people with impaired hand function.



#### Adapts to Your Needs

A control pad on the wrist lets users switch between three customizable profiles or activate the glove for grip assistance. Profiles can be tailored for different tasks. The profiles are configured via an app.

## CanDo® Inflatable Exercise Balls-Standard

Inflatable Exercise Balls

Origin: USA



- Use to improve balance, coordination, flexibility and strength, or just for fun
- Used by therapists as an aid for vestibular movement and equilibrium therapy, by fitness professionals as an aid for stretching and strength programs, as seats, and even by pregnant women in birthing classes
- Non-slip surface is ribbed for extra security

- Underinflate to give a soft, mushy feel or fully inflate for a firm, bouncy feel
- Supports up to 300 lb (136kg)

Retail Box		
Cm	Color	Code
45	Yellow	30-1801B
55	Orange	30-1802B
65	Green	30-1803B
75	Red	30-1804B
85	Blue	30-1805B

Polybagged		
Cm	Color	Code
30	Blue	30-1800
45	Yellow	30-1801
55	Orange	30-1802
65	Green	30-1803
75	Red	30-1804
85	Blue	30-1805
95	Red	30-1806
105	Blue	30-1841
120	Orange	30-1807
150	Green	30-1808

## CanDo® Closed Cell Exercise Mats

Pediatric Mats

Origin: USA  
Certificate: CE

- Cushioned exercise mats offer superior comfort, quality and durability, yet are lightweight and sturdy enough for commercial use in schools and childrens' facilities
- Colorful, soft, resilient closed-cell foam
- Waterproof and sanitary; can even be used in pool
- Dimensions (L x W x H) 72" x 24" x 0.6"

Colour	Each
Yellow	30-2310Y





## SUP-R BAND

Latex Free Exercise Band

Origin: USA

### Features

- Choose the Twin-Pak® and save
- Use with all dispenser systems
- Textured surface is easy-to-grip
- Powderless band leaves no mess

### Twin-Pak 100 Yard Dispenser Box (2) 50-Yard Boxes

Tan	xx-light	10-6330
Yellow	x-light	10-6331
Red	light	10-6332
Green	medium	10-6333
Blue	heavy	10-6334
Black	x-heavy	10-6335
Silver	xx-heavy	10-6336
Gold	xxx-heavy	10-6337
Yel Red Grn Blu Blk	set (1 ea)	10-6338



### 50 Yard Dispenser Box

Tan	xx-light	10-6320
Yellow	x-light	10-6321
Red	light	10-6322
Green	medium	10-6323
Blue	heavy	10-6324
Black	x-heavy	10-6325
Silver	xx-heavy	10-6326
Gold	xxx-heavy	10-6327
Yel Red Grn Blu Blk	set (1 ea)	10-6328
Set of 5 w/ Dispens-a-Band rack		10-6329

### 6-Yard Dispenser Box

Tan	xx-light	10-6310
Yellow	x-light	10-6311
Red	light	10-6312
Green	medium	10-6313
Blue	heavy	10-6314
Black	x-heavy	10-6315
Silver	xx-heavy	10-6316
Gold	xxx-heavy	10-6317
Yel Red Grn Blu Blk	set (1 ea)	10-6318

challenge level	tan 0	yellow 1	red 2	green 3	blue 4	black 5	silver 6	gold 7
	easy			moderate				difficult



Set With Dispens-A-Band® Rack

## CanDo® Multi - Grip Exerciser

Easy Grip Loops Allow For Unlimited Exercise Options

Origin: USA

### Features

- Each exerciser contains both large (6") and small (2") grips. Use small grips as anchor or for finger exercises
- Loop around hands and feet to perform high impact upper and lower body exercises
- Cloth exerciser contains no latex, no scent, and no powder coating
- Can be washed and dried
- Easily change resistance by moving 1 loop up or down the exerciser, or using a different resistance
- Progressively exercise with different resistance levels
- Available in 9-loop (6 ft. exercisers) and 135-loop 15 exercisers (30 yd. rolls)



30 yard roll  
15 exercisers



9 loop (6 foot)  
exerciser

		9 Loops, 6 Foot Exerciser	9 Loops, 6 Foot Exerciser (24 Ea)	135 Loops, 15 Exercisers, 30 Yard Roll
Tan	xx-light	10-6600	10-6600-24	10-6610
Yellow	x-light	10-6601	10-6601-24	10-6611
Red	light	10-6602	10-6602-24	10-6612
Green	medium	10-6603	10-6603-24	10-6613
Blue	heavy	10-6604	10-6604-24	10-6614
Black	x-heavy	10-6605	10-6605-24	10-6615
Silver	xx-heavy	10-6606	10-6606-24	10-6616
Gold	xxx-heavy	10-6607	10-6607-24	10-6617
Yel Red Grn Blu Blk	5pc set	10-6608	10-6608-24	10-6618

## CanDo® Exercise Band Loops

High Quality Latex Exercise Band Pre-Formed Into A Loop

Origin: USA

### Features

- Available in 3 lengths to perform upper, lower and full body exercises
- Lengths represent loop when flat; loops are 3" wide
- Full Body Sets contain 10", 15" and 30" loop for total body workout
- Black, silver and gold only offered in CanDo®

#### Bulk Packaging (10 Pack)

		10" Loop	15" Loop	30" Loop
Yellow	x-light	10-5251-10	10-5261-10	10-5291-10
Red	light	10-5252-10	10-5262-10	10-5292-10
Green	medium	10-5253-10	10-5263-10	10-5293-10
Blue	heavy	10-5254-10	10-5264-10	10-5294-10
Black	x-heavy	10-5255-10	10-5265-10	10-5295-10
Silver	xx-heavy	10-5256-10	10-5266-10	10-5296-10
Gold	xxx-heavy	10-5257-10	10-5267-10	10-5297-10
Yel Red Grn Blu Blk	set (1 ea)	10-5259-10	10-5269-10	10-5299-10



30" loop 0" loop  
(full body) (11 body)

## CanDo® Inflatable Exercise Rolls

Inflatable Exercise Rolls &amp; Balls

Origin: USA

### 1 Saddle Roll

- Easier to control because motion is limited to one direction
- Aids in developing muscles, coordination and balance
- Latex and phthalate-free
- Supports up to 300 lb (136 kg)

Cm-Dia X L	In-Dia x L	Color	Code
40 x 90	16 x 35	Yellow	30-1725
50 x 100	20 x 39	Orange	30-1726
60 x 110	24 x 43	Green	30-1727
70 x 120	28 x 47	Red	30-1728
80 x 130	32 x 51	Blue	30-1729



### 2 Sensi-Saddle Roll

Cm-Dia X L	In-Dia x L	Color	Code
40 x 90	16 x 35	Yellow	30-1735
50 x 100	20 x 39	Orange	30-1736
60 x 110	24 x 43	Green	30-1737



### 3 Straight Roll

Cm-Dia X L	In-Dia x L	Color	Code
40 x 90	16 x 35	Yellow	30-1780
50 x 100	20 x 43	Orange	30-1781
60 x 110	24 x 53	Green	30-1782

## CanDo® Inflatable Balance Discs

Vestibular Balance Discs &amp; Wedges

Origin: USA

- When in a seated position, inflatable balance disc mimics the movement and shape of an inflatable ball
- Use on floor as a standing disc for balance training, proprioception and strengthening of the core muscles
- Use with children who have trouble sitting still to increase attention and focus
- Disc has one flat side and one nubby side for tactile feedback and stimulation

	24" / 60cm Diameter	14" / 35cm Diameter
Black	30-1868BLK	30-1870BLK
Blue	30-1868B	30-1870B
Green	30-1868G	30-1870G
Red	30-1868R	30-1870R
Yellow	30-1868Y	30-1870Y
Silver	—	30-1870S
Cover (Not Shown)	30-1867	—



## CanDo® Web Hand Exercisers

Hand & Wrist Exercisers

Origin: USA

- Perform flexion, extension, opposition and supination exercises
- Modify exercises by adjusting the hand position or depth of finger insertion
- Modify resistance by switching to a different resistance color web
- Available in both 14" and 7" diameter size; latex and latex-free
- Resistance levels are color-coded



Color		Large 14" Diameter		Small 7" Diameter	
		Low Powder	Latex-Free	Low Powder	Latex-Free
Tan	xx-light	10-0850	10-0870	10-0860	10-0880
Yellow	x-light	10-0851	10-0871	10-0861	10-0881
Red	light	10-0852	10-0872	10-0862	10-0882
Green	medium	10-0853	10-0873	10-0863	10-0883
Blue	heavy	10-0854	10-0874	10-0864	10-0884
Black	x-heavy	10-0855	10-0875	10-0865	10-0885
Tan Yel Red Grn Blu Blk Set		10-0856	10-0876	10-0866	10-0886
Multi-Resistance 14" Large Diameter, Latex					
Yellow Green	x-light/medium	10-0857			
Red Blue	light/heavy	10-0859			

## CanDo® Twist-n-Bend Exercisers

Hand & Wrist Exercisers

Origin: USA

- Change resistance by varying distance between hands and/or by using a different resistance bar
- TWIST:** with both hands on the bar, twist by flexing one wrist and extending the other
- BEND:** place both hands on or near the ends of the bar; keep your hands in position and exert force to bend the bar
- TWIST-N-BEND:** combine Twist and Bend methods
- SHAKE:** place hand(s) at center or end, shake the bar



Color		12" Long
Tan	xx-light	71-0120
Yellow	x-light	71-0121
Red	light	71-0122
Green	medium	71-0123
Blue	heavy	71-0124
Black	x-heavy	71-0125
Tan Yel Red Grn Blu Blk		30 pieces set, 5 of each colour 71-0110

## CanDo® Digi-Extend Hand Exercisers

Hand Exercisers

Origin: USA

- Facilitates extension and flexion of all fingers
- Develop isolated finger strength, flexibility and coordination
- 4 color-coded resistance band levels included (tan, yellow, red, green) allow each muscle or joint to work at a comfortable level and build as increased strength is achieved
- For rehabilitation of carpal tunnel, arthritis, stroke, fractures, tendon injury, nerve lacerations, tennis elbow and more



Exerciser	
10-0775	CanDo Digi-Extend Hand Exerciser with 14 Bands (5 Tan, 4 Yellow, 3 Red, 2 Green)
Clinic Set	
10-0757	CanDo Digi-Extend Exercisers, Set of 4, 56 Bands Total, Metal Stand
Metal Stand	
10-0756	Accessory, 4 Position Metal Stand Only

Additional Latex-Free Bands			
Replacement Resistance Bands - 25 Each			
10-1850	Tan	xx-light	
10-1851	Yellow	x-light	
10-1852	Red	light	
10-1853	Green	medium	
10-1854	Blue	heavy	
25 Latex-Free Bands - 5 Each Color			
10-0756	Stand Only	All Colors	set

## CanDo® Digi-Flex Hand Exercisers

Hand &amp; Finger Exercisers In 8 Strength Levels

Origin: USA

- Upper and lower spring design allows more flexibility for all hand motions
- Develops isolated finger strength, flexibility and coordination
- For rehabilitation of carpal tunnel, arthritis, stroke, fractures, tendon injury, nerve lacerations, tennis elbow and more

Color	Force To Compress Finger/Hand (Lb)	Code
Tan	xx-light	10-0751
Yellow	x-light	10-0740
Red	light	10-0741
Green	medium	10-0742
Blue	heavy	10-0743
Black	x-heavy	10-0744
Silver	xx-heavy	10-0752
Gold	xxx-heavy	10-0753
Yel Red Grn Blu Blk	set of 5 (1 ea)	10-0745
Yel Red Grn Blu Blk	set with plastic rack	10-0746
Yel Red Grn Blu Blk	set with metal rack	10-0749
Tan Yel Red Grn Blu Blk Sil Gld	set of 8 (1 ea)	10-0758
Tan Yel Red Grn Blu Blk Sil Gld	set of 8 with metal rack	10-0759





## CanDo® TheraPutty® Standard Exercise Putty

TheraPutty® Hand Exercise Material

Origin: USA

- The standard in resistive hand exercise material
- Each color-coded putty has a different firmness ranging from xx-soft to x-firm
- Available in the Preferred® color sequence
- Available in convenient resealable plastic containers (2, 3, 4 and 6 ounce) for patients to use at home and a choice of larger containers (1, 5 and 50 pound) for the hospital or clinic
- CanDo® TheraPutty® material is gluten, casein and latex-free; non-toxic if swallowed
- CanDo® TheraPutty® material putty does not dry out



	Standard Putty - Each						Putty - Set	
	Tan	Yellow	Red	Green	Blue	Black	Complete Anti-Microbial Set	
	xx-soft	x-soft	soft	medium	firm	x-firm	1 each level (6 total) tan yel red grn blu blk	
2 oz	10-0956	10-0900	10-0901	10-0902	10-0903	10-1466	2 oz	10-1480
3 oz	10-0957	10-0967	10-0968	10-0969	10-0970	10-1467	3 oz	10-1481
4 oz	10-0958	10-0905	10-0906	10-0907	10-0908	10-1468	4 oz	10-1482
6 oz	10-0959	10-0910	10-0911	10-0912	10-0913	10-1469	6 oz	10-1483
1 lb	10-0994	10-0918	10-0919	10-0920	10-0921	10-1470	1 lb	10-1484
5 lb	10-0995	10-0923	10-0924	10-0925	10-0926	10-1471	5 lb	10-1485
50 lb	10-1460	10-1461	10-1462	10-1463	10-1464	10-1473		

## CanDo® Theraputty Antimicrobial Exercise Putty

Origin: USA

- Same resistance levels and sizes as standard CanDo® TheraPutty® hand exercise material
- Available in the Preferred® color sequence
- CanDo® TheraPuttyPlus® anti-microbial material is gluten, casein and latex-free
- Non-toxic anti-microbials are used to kill organisms

	Standard Putty - Each						Putty - Set	
	Tan	Yellow	Red	Green	Blue	Black	Complete Anti-Microbial Set	
	xx-soft	x-soft	soft	medium	firm	x-firm	1 each level (6 total) tan yel red grn blu blk	
2 oz	10-2600	10-2601	10-2602	10-2603	10-2604	10-2605	2 oz	10-2606
3 oz	10-2610	10-2611	10-2612	10-2613	10-2614	10-2615	3 oz	10-2616
4 oz	10-2620	10-2621	10-2622	10-2623	10-2624	10-2625	4 oz	10-2626
6 oz	10-2630	10-2631	10-2632	10-2633	10-2634	10-2635	6 oz	10-2636
1 lb	10-2640	10-2641	10-2642	10-2643	10-2644	10-2645	1 lb	10-2646
5 lb	10-2650	10-2651	10-2652	10-2653	10-2654	10-2655	5 lb	10-2656
50 lb	10-2660	10-2661	10-2662	10-2663	10-2664	10-2665		



**Disclaimer:**

For binding information on Hocoma/Motek/DIH products, please refer to [www.hocoma.com](http://www.hocoma.com) / [www.motekmedical.com](http://www.motekmedical.com) / [www.dih.com](http://www.dih.com) . Hocoma/Motek/DIH does not take any responsibility for the accuracy, completeness and quality of the information provided on this catalog. Liability claims regarding damage caused by the use of any information provided, including any kind of information which is incomplete or incorrect, will therefore be rejected.